

DETECTO®

495 SERIES

STATIONARY WHEELCHAIR SCALE



Long-Lasting Weighing Accuracy Without Any Power Required

- 450 lb Capacity
- 4 oz Accuracy
- Beam is Readable From Both Sides
- Die-cast, Easy-to-read Beam
- Skid-proof Black Mat
- Spacious Weighing Platform
- Ramp Included
- Sturdy Steel Construction

DETECTO SCALES MAKE WEIGHING WHEELCHAIR-BOUND PATIENTS EASIER AND FASTER

Eye-Level Weight Display for Healthcare Professionals



TIME-TESTED WEIGHBEAM

DETECTO's mechanical weighbeam is world-renowned for its high accuracy and long-term dependability over decades of use.

READABLE FROM BOTH SIDES

The weighbeam is imprinted on both sides, so the nurse can read the patient's weight from either side, front or back.

- No Power Required: Mechanical Weighbeam Style
- Detachable Ramp
- Long-Lasting, Clinical-White Powder Paint Finish



This heavy-duty stationary wheelchair scale features a detachable ramp for smooth entry and exit. The 495 features a padded handrail for ambulatory patients and DETECTO's world-renowned die-cast, dual-sided, eye-level weighbeam graduated in 4 oz increments. The weigh beam is readable from either side.

Input from customers specializing in the care of ambulatory patients provided information which was transferred into features and benefits of this scale, evident the first time you use it! The weighbeam is located at a height above the platform and wheelchair location where a healthcare professional can perform the weighing operation easily and comfortably.

Model	495
Capacity	450 lb x 4 oz
Platform Size	30 in W x 26 in D x 2 in H / 76 cm W x 66 cm D x 4 cm H
Shipping Weight	187 lb / 85 kg
Beam Type	Mechanical
Design	Equipped with wheelchair ramp



DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.



SOLD BY:

102 East Daugherty St., Webb City, MO 64870 USA
 Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com