

SEATED, STANDING,



The 6560 and 6570's dual-access ramps go all the way to the ground and allow maximum convenience for loading and unloading patients from either side of the scale. These unique ramps require only half the effort for a healthcare professional to roll a wheelchair-bound patient onto the scale than competitor's models.

- Clinical-grade accuracy to 0.2 lb / 0.1 kg
- Quick and easy mobility for transport or storage
- Up to 99 tares (optional) to store wheelchair weights for efficient patient weighing
- Hold/Release key to hold weight after leaving the scale for better patient care
- Interface with a PC or printer via serial, USB or optional Wi-Fi or wired Ethernet for EMR
- Patient ID barcode scanning sends 14-digit patient ID to indicator for EMR/EHR records

SE

DIGITAL FILTERING

elchA

StableSEN

Model

6570



Handrail for Patient Stability The 6560 and

Curved

6570's elegant curved handrail provides patient support while weighing standing on the platform.

Single-Pass Weighing Saves Time

Enjoy the time-saving benefits of these scales by using the convenient pushbutton or keypad tare to remove the wheelchair weight.



Fold-up **Chair Seat**

The 6570's chair seat easily folds up when not in use, allowing wheelchair weighing (not available on 6560).

The scale's columns fold down and latch securely for easy transport and vertical foldup storage requiring minimal floor space.

Easily-Collapsible Column

By sliding the release lever, you can easily fold down the columns for easy transport or vertical storage.

RULTIPLE CLINICAL DEPARTMENTS

Digital Multifunctional Medical Scales with Time-Saving and Patient Comfort Benefits

DETECTO's 6560 and 6570 portable wheelchair scales offer time-saving and patient comfort benefits. The scales features 1,000 lb / 450 kg capacity, optional Wi-Fi for EMR/EHR, vertical foldup storage, two-way built-in ramps for accessibility from both sides, large 32 in x 36 in / 81 cm x 91 cm platform, clinical-grade accuracy to 0.2 lb / 0.1 kg, battery or AC power, and up to 99 tares (optional) to store wheelchair weights for efficient patient weighing. Enjoy the time-saving benefits of these scales by using the convenient pushbutton or keypad tare to remove the wheelchair weight. OIML EC class III models are available (see back page).

EFFICIENT MEASUREMENT



The 6560 and 6570's MedVue® weight indicator is the clinical industry's most advanced, featuring EMR/EHR-ready serial, USB, and optional Wi-Fi connectivity, BMI calculation, yellow "quick keys" for basic weighing, and a multiline LCD display. The scale is powered by 6 C batteries (not included) or optional AC adapter.

Low-Profile, Thin Platform

The platform is only 2.2 inches/56 mm high for easy patient accessibility and comfort while weighing either standing or in a wheelchair.



The scale's columns fold down and latch

securely for easy transport and vertical fold-up storage requiring minimal

floor space.

	Capacity	6560 1,000 lb x 0.2 lb /	6560KGEU 300 kg x 0.1 kg or	6570 1,000 lb x 0.2 lb /	6570KGEU 300 kg x 0.1 kg or
ready REGISTERED 9001	OIML EC Class III Approved	450 kg x 0.1 kg	350 kg x 0.1 kg 01	450 kg x 0.1 kg	350 kg x 0.2 kg
	Platform Size	32 in x 36 in x 2.2 in / 81 cm x 91 cm x 5.6 cm			
	Keypad	19 key, color-coded, sealed membrane			
	Connectivity	RS232 serial, USB and optional wired and Wi-Fi Ethernet			
	LCD Display Dimensions	3.5 in W x 3 in H / 89 mm W x 76 mm H backlit blue LCD			
÷13485	Power	Optional model MV1PWR AC adapter or 6 C batteries (not included)			
	Resolution	5,000 divisions, max.			
	Operating Temperature	14° to 104°F (-10° to 40°C)			
	Ramp Size		36 in (91.4 cm) W x 7.9	in (20.2 cm) L eacl	n
	Display	5 digit, 7 segment, 7/8 inch (22.23 mm) high LCD weight digits			
	Country of Origin	•••••••••••••••••••••••••••••••••••••••			
	UPC Code	809161302202	809161304503	809161302301	809161304602
	Net Weight	116 lb / 53 kg	116 lb / 53 kg	140 lb / 63 kg	140 lb / 63 kg
	Shipping Weight	126 lb / 57 kg	126 lb / 57 kg	150 lb / 68 kg	150 lb / 68 kg

DETECTO reserves the right to improve, enhance,

SOLD BY:

or modify features and specifications without prior notice.



203 E. Daugberty, Webb City, MO 64870 USA Pb: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153

Mobile: m.detecto.com