

# 6745

## DIGITAL BABY SCALE



- Capacity: 30 lb x 0.1 oz / 15 kg x 0.005 kg
- Lock-down weighing tray can be removed for easy cleaning
- Zero button allows you to remove blanket tare weights
- Metric lockout available in standard model and also kilograms-only model
- Baby scale carts also available from DETECTO

# WEIGH THE SMALLEST OF PATIENTS

WITH HIGH ACCURACY AND DEPENDABILITY



## Designed for Measuring Active Babies

- Lock-in weight feature which compensates for typical newborn movement
- Spacious weighing tray: 26 in W x 12 in D / 66 cm W x 31 cm D
- 6 C batteries (not included) or AC power (AC adapter included)
- Large, easy-to-read 1-in / 25-mm high LCD readouts
- Integral measuring tape included allowing weight and height measurements



DETECTO's 6745 baby scale was designed for use in the NICU, nursery, delivery room and pediatric practices where high accuracy and dependability is paramount for weighing the smallest of patients. Step-by-step instructions are displayed right on the keypad for ease of use.

If you wish to place a blanket or pad on the baby tray, you can press the ZERO key to reset the weight display to zero. This allows you to remove the blanket tare weight. Toggle easily between pounds and kilograms by pressing the UNITS key (metric lockout is also available in set-up).



## PC OR PRINTER CONNECTIVITY FOR PATIENT RECORDS

The 6745's RS232 port allows it to be connected to a printer, such as the DETECTO P50 or P150 which include the cable needed to communicate with the scale, or with a PC for outputting to EMR/EHR software for patient medical records.

# DIGITAL BABY SCALE WITH FINE .1 OZ .005 KG ACCURACY



The locking pin on the front of the weighing tray allows you to quickly release it for easy cleaning and then lock back into place.

## From Neonatal to Pediatric Use

The ergonomically-shaped tray with molded sides makes it easy to weigh larger babies and keeps them safe and secure while measuring.

The DETECTO 6745 baby scale offers an integral bubble level and four adjustable feet which allow you to quickly find a level weighing surface. The 6745 can use six C size Alkaline, Ni-Cad or NiMH batteries (not included) or the universal medical-grade AC adapter (included). The scale can operate for 200 hours of continuous use when using alkaline batteries or 50 hours of continuous use with fully charged Ni-Cad or NiMH batteries.

With strict dosage measurement requirements needed for hospitals, the 6745 offers a variety of weighing units settings:

- Pounds and Ounces Only
- Kilograms Only
- Pounds and Ounces / Kilograms
- Kilograms / Pounds and Ounces
- Or you may purchase the kilograms-only model 6745KG

- Built-in battery charger
- Color-coded, easy-to-use keypad
- Serial output allows the scale to be connected to a printer or PC
- 200 hours of continuous use with Alkaline batteries
- Sleep and Auto Shutoff modes extend the battery life
- Adjustable lock threshold
- Acceleration of gravity compensation

## BABY SCALE CARTS

Rolling baby scale carts with locking casters are also available from DETECTO to easily move the 6745 from location to location.



Model SPBT1728 stainless steel cart with open shelf.



Model RC33669WHT steel cart with powder coat paint and 5 keyed-lock drawers.



# 6745

## DIGITAL BABY SCALE



### CONNECTIVITY PORT

The 6745 baby scale features RS232 serial output for connectivity with a PC or printer for patient records.



### FLEXIBLE POWER

Your choice of rechargeable NiCad batteries, standard Alkaline batteries, or UL-listed AC adapter.

MODEL	6745	6745KG
Capacity	30 lb x 0.1 oz / 15 kg x 0.005 kg	15 kg x 0.005 kg
Tray Size	26 in W x 12.2 in D / 66 cm W x 31 cm D	
Overall Size	26 in W x 14.2 in D x 7.5 in H / 66 cm W x 36 cm D x 19 cm H	
Display	6-digit, 7-segment 1.0 in / 25 mm high transflective LCD	
Power	6 C size Alkaline, Ni-Cad or NiMH batteries (not included) or 100-240VAC, 12VDC/1.2A medical UL-approved AC adapter (included)*	
Weight Units	Pounds and Ounces / Kilograms	Kilograms
Integrated Measuring Tape	0 - 22 in / 0 - 56 cm	
Connectivity Port	RS232 serial	
Keys	ON, ZERO, UNITS, PRINT, OFF	
Net Weight	18 lb / 8 kg	
Shipping Weight	22 lb / 10 kg	
Country of Origin	USA	
UPC Code	809161140705	809161302806

\* The scale can operate for 200 hours of continuous use when using alkaline batteries, or 50 hours of continuous use with fully charged Ni-Cad or NiMH batteries.

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.



203 E. Daugherty, Webb City, MO 64870 USA  
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153

Mobile: [m.detecto.com](http://m.detecto.com)

### SOLD BY: