

ñn

DIGITAL BABY SCALES WITH INCLINED SHELL SEAT

Infants readily accept the inclined seating of DETECTO's 8432-CH series digital pediatric scales. The built-in, shell-shaped seat is made of plastic material and fun for children to sit in while weighing. With pound and kilogram readings and heavyduty steel base and understructure for support and stability, these pediatric scales have a lot to offer.

- Digital speed and accuracy
- 44 lb / 20 kg capacity
- Clinical accuracy to .05 lb / .02 kg
- Simple-to-use weight indicator
- Battery or AC power
- Solid steel base for steady assurance while weighing
- Pushbutton tare removes blanket weights
- Units locking into pounds or kilograms for safe dosages

| | | BY SC D SHEL | | The DETECTO 8432-CH pediatric scale with in seat is ideal for use in l clinics and medical pra weight measurement infants and tots. The sc an easy-to-clean, seashell-shaped plastic seat for optimal comfort of the child and a durable, sturdy | clined chair health actices for of |
|--------------------|---|--------------------------------------|------------------|---|---|
| | | | | enameled steel scale base. | 11 |
| | | | | | |
| | | | | | |
| | | | | | |
| Model | 8432-CH | 8432-CH-AC | 8432KG-CH | | |
| Capacity | 44 lb x 0.05 lb / 20 kg x 0.02 kg | 44 lb x 0.05 lb / 20 kg x 0.02 kg | 20 kg x 0.02 kg | | / |
| AC Adapter | | • | | | |
| Scale Height | 26 in / 66 cm | | 15. 75 in | - | |
| Display | Six-digit, seven segment, 0.7 inch (17.8 mm) high LCD | | '.8 mm) high LCD | | (|
| Power | Six "AA" Alkaline, Ni-Cad or NiMH batteries (not included) OR an optional 100 to 240V AC 50/60Hz 12V DC 1A wall plug-in UL/CSA listed AC power adapter (Cardinal part number 6800-1045). Models are available with the AC adapter included. | | | | |
| Keyboard | Membrane type with 5 keys and 4 directional arrows | | | | |
| Connectivity | RS232 serial port | | | | |
| Shipping Weight | 33 lb / 15 kg | | | 11 in 28 cm 12.25 in 31 cm | |



The 8432-CH has a lock-in weight feature which compensates for typical baby movement. The crisp LCD features 0.75 inch high digits for easy readouts.



DETECTO's plastic shellshaped inclined seat on the 8432-CH series baby scales makes weighing fun for pediatric patients. The weight indicator is conveniently situated over the top of the shell seat, so it's easy to read the weight

while the baby is seated.



102 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153 www.Detecto.com DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

SOLD BY: