

WET DIAPER SCALE

AP-4KD



- Easy-to-read, 0.5 in red digital LED display
- Compact size fits any space
- Removable stainless steel tray for easy cleaning
- Pushbutton tare

Model AP-4KD WET DIAPER SCALE

Accurate, compact, and easy-to-use... DETECTO's wet diaper scale has been expertly engineered to help you make the most of your valuable time.

ONE COMPACT SCALE



REMARKABLE ACCURACY QUALITY EASE OF USE

ACCURACY

This compact scale features a convenient 4,000 gram capacity and is accurate to 1 gram.

QUALITY


DETECTO is the largest medical scale manufacturer in the United States providing clinical-grade accuracy and dependability for over a century.

EASE OF USE

DETECTO's model AP-4KD has been designed to meet your every need in measuring fluid output, whether from a diaper or sanitary pad. The stainless steel construction and removable tray makes cleanup a snap... a major consideration in fast-paced hygienic work environments. An instantaneous 0.5" high red digit display is easy to read when time is of the essence. Handy pushbutton tare simplifies weight determinations.

Capacity	4,000 g x 1 g
Dimensions	7.0 in L. x 7.5 in W. x 2.8 in H. (17.8 cm x 19.1 cm x 7.1 cm)
Removable Tray	9.0 in x 6.0 in stainless steel (22.9 cm x 15.2 cm)
Shipping Weight	5.0 lb / 2.3 kg
Power	115 VAC, 50/60 Hz or 230 VAC, 50/60 Hz; AC adapter with 6 ft power cord included
Tare	Pushbutton tare

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

 **DETECTO**[®]
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Mobile: m.detecto.com

SOLD BY: