

# ALL-DIGITAL PORTABLE WEIGHT AND MEASUREMENT



## Sonaris Touchless Sonar Stadiometer

- Measures 1 to 9 ft (30 to 275 cm)
- Pinpoint accuracy to 0.2 in and 5 mm graduations



DETECTO  
Model SONARIS

- Instant digital height readouts without pushing a single button
- May be wall mounted or set on a countertop

- Arrives fully assembled and ready for measurement
- Auto calibrate (or manual calibration) to current position for completely portable use
- Digitally measures small children up to tall adults



HOSPITALS  
SCHOOL HEALTH  
MEDICAL EXAMS  
PATIENTS' HOMES • DOCTORS' OFFICES




## APEX-RI Remote Indicator Portable Scale

- 600 lb x 0.2 lb / 300 kg x 0.1 kg
- Spacious 17 x 17 in / 43 x 43 cm platform



- Portable with built-in carrying handle
- Textured platform surface for safety
- Remote indicator for versatile display placement
- Body Mass Index
- 12VDC AC power adapter included on -AC models



DETECTO  
Model APEX-RI

**DETECTO**