

AMBULATORY SURGERY CENTER PRODUCTS

Scales and related products for ASCs
Modern innovative scales for modern healthcare facilities



8550

6880

WPS12UT

MB130 with DLM

7550

solace
IN-FLOOR DIALYSIS SCALES

ID-7248SH-855RMP

RESCUE
medical carts

apex
CLINICAL SCALES

Stainless Steel
Step-On Cans

RC33669BLU-L
RC33669RED-L

6570

APEX-SH

WASTE MATE
Waste Mate Plastic
Waste Receptacles

DETECTO

ACCURATE AND RELIABLE SCALES AND EQUIPMENT FOR OUTPATIENT OPERATING ROOMS

Offering the ultimate in patient satisfaction and accuracy

DETECTO manufactures the widest variety of medical weighing products on the market for every aspect of ambulatory surgery center patient care... from seated, standing, in-bed, pediatric, and wheelchair patient weighing, every type of scale is covered in DETECTO's clinical product portfolio.

7550

Wall-Mount Fold-Up Scale



The 7550 wall-mount wheelchair scale features a fold-up design that saves space when not in use. It can be installed in high-traffic areas for convenient weighing, and fold up when finished weighing. The scale accommodates most wheelchairs, and yet requires very little floor space when folded up.

6570

Portable Wheelchair Scale with Handrail and Seat



The 3-in-1 6570 scale provides seated, standing and wheelchair weighing. The fold-away seat and stable handrail offer maximum flexibility without sacrificing portability. The large 32 in x 36 in platform has two-way built-in ramps for maximum accessibility and 1,000 lb/450 kg capacity.

ID-7248SH-855RMP

Solace In-Floor Flush-Mounted Dialysis Scale



Solace series in-floor dialysis scales are flush-mounted platform scales that feature stainless steel decks and include a pit frame, trim ring, and interface cable to the included multi-color, touchscreen LCD indicator and printer. The flush surface eliminates tripping hazards and accessibility issues.

APEX-SH

Digital Physician Scale with Touchless Sonar Height Rod



The apex® digital scale with sonar height rod lets you measure height and weight without pushing a single button. The large .75-inch-high LCD readout is easy to view. The scale provides BMI, runs on AC or battery power, and comes with a variety of connectivity options including Wi-Fi/Bluetooth.

WPS12UT

Lap Sponge Scale is Submersible Underwater



The WPS12UT digital scale with utility bowl is IP67-rated for water protection and can be completely submersed underwater for easy cleaning. This stainless steel scale is ideal for weighing soiled or wet diapers, chux pads, lap sponges, and other materials. The 4.6 in / 11.7 cm deep bowl aids sample retention.

MB130 with DLM

Digital Pediatric Scale with Digital Length Measuring



This flexible scale lets you weigh patients from neonatal up to toddler with its 44lb/20kg capacity. The extra-large weighing tray keeps babies safe and comfortable, and the 5 gram accuracy is great for weighing your smallest patients. A measuring tape is built in, with a digital length measuring rod available.

8550

Stretcher Scales with Column-Mounted Indicator



DETECTO's Stretcher Scale is perfect for quickly weighing emergency arrivals, with a portable design that makes it easy for anyone to move. The column display is easy to read, with a Hold/Release key to hold the patient's weight after leaving the scale for better patient care.

6880

Chair Scale with Liftaway Armrests



DETECTO's 6880 chair scale is the simple solution for patients with decreased mobility and features 550 lb / 250 kg capacity, dual wheels for transport, smartphone-style indicator, liftaway padded armrests and footrests, and battery or AC power. Wi-Fi and Bluetooth wireless connectivity models are available.

RC33669BLU-L RC33669RED-L

RESCUE Steel Medical Carts with Customizable Accessories



DETECTO's Rescue medical carts are sturdy and reliable... built to make sure you have the supplies you need where and when you need them. Carts have your choice of five or six drawers that can be secured with a key or EMG breakaway lock and the ABS countertop is easy to clean.

C-48

Step-On Cans, the Rugged Waste Receptacle



Step-On Cans are perfect for ASCs and come equipped with a heavy-duty foot pedal and linkage and a bag securing mechanism, they come in a variety of five different sizes and capacities and two different materials: stainless steel or a baked epoxy finish (available in red or white).

WM-18RED

Waste Mate Plastic Waste Receptacles



Available in four sizes and red and beige colors, DETECTO's Waste Mate commercial-grade plastic waste receptacles feature an easy-to-use foot pedal that flips open the lid for no-touch trash disposal. The integrated internal clamps secure the trash bag to the rim of the can to prevent it from slipping into the can.

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

SOLD BY:



102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com