

DETECTO®

Model
SPBT-1728

Rolling Baby Scale Cart



Detecto's stainless steel rolling baby scale cart allows the weighing to come to the patient.

The smooth gliding 2.5" diameter rubolene wheels and 27.75"W x 16.75"D x 33.4" H cart dimensions allow it to easily fit through doorways into patient rooms. A lower shelf with three height adjustment levels allows clinical staff to store patient records and instruments onboard the cart. The spacious shelf space is the perfect size for baby scales, such as Detecto's model 6745 (pictured), which is battery operated for complete portability.