

BARIATRIC SCALES

1028

DETE

CTO

featuring DIGITAL Height Rods

DIGITAL SPEED AND ACCURACY

Heavy-duty, yet lightweight, the 6857DHR bariatric scale's balanced design and integral wheels make it easy to move. The 6857DHR is engineered specifically for weighing obese and unsteady patients and features a roomy, spacious 2-ft-square platform.

DETECTO Model 6857DHR

Bulletin No. C227B

DETECTO DIGITAL BARIATRIC SCALES -- WITH DIGITAL HEIGHT RODS --



DIGITAL WEIGHT • HEIGHT • BMI

The 6857DHR and 6854DHR bariatric scales are designed for optimum performance and display weight, height, and BMI simultaneously (shown below) by

VERSATILE CONNECTIVITY

Detecto's state-of-the-art MedVue[®] digital weight indicator includes DB9 serial, USB, and optional wired or Wi-Fi Ethernet output for data transfer to a printer or PC-based EMR (electronic medical record) software. With the time/date feature, the scale can send the weight, time, date, BMI (Body Mass Index), and height directly to a printer or a PC for record keeping. The 350-degree swivel mounting bracket allows it to be positioned wherever it's easiest for viewing. simply raising the digital height rod while a patient is being weighed on the scale. There are no keys to push, no operational sequences to memorize, it's fully automatic. These bariatric scales are EMRready with USB, serial, and optional Wi-Fi. The ability to print or save weight data with these scales is a real time-saving device for clinics, doctor's offices, long-term care facilities, and other medical institutions.



Digital measurement speed and accuracy
Stainless steel handrails for assurance

Integral wheels for quick and easy transport

assurance 💫 🕨 Height rod measures up to 6'7" / 200 cm

	DETECTO DIGITAL BARIATRIC SCALE SPECIFICATIONS			
DETECTO Model 6857DHR		6857DHR	6854DHR	
	CAPACITIES:	1,000 lb x 0.2 lb / 450 kg x 0.1 kg	600 lb x 0.2 lb / 270 kg x 0.1 kg	- dut-word
	PLATFORM SIZE:	24" x 24" / 61 cm x 61 cm	18" x 14" / 46 cm x 36 cm	
	SHIPPING WEIGHT:	85 lb / 39 kg	74 lb / 34 kg	
	DIGITAL HEIGHT ROD MEASUREMENTS:	43" / 110 cm - 79" / 200 cm	same	
	POWER:	6 AA size alkaline batteries (not included) or optional AC adapter	same	
	DISPLAY:	Liquid Crystal, 0.875"/22 mm high weight readouts	same	
	CONNECTIVITY:	Serial, USB, and optional Wi-Fi	same	DETECTO Model 6854DHF

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.



Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001 www.Detecto.com Sold By: