

# CHAIR SCALES

Seated healthcare scales for weighing patients that are frail, handicapped, on dialysis, or having mobility issues



475

6880

6475

**DETECTO**<sup>®</sup>

# DIGITAL AND MECHANICAL PORTABLE CHAIR SCALES WITH **LIFT-AWAY** ARMRESTS AND FOOTRESTS

DETECTO's chair scales are a welcome addition to any healthcare facility when weighing patients that are frail, handicapped, on dialysis, or having mobility issues. DETECTO's convenient chair scales can save a great deal of time and effort when a healthcare facility does not have the proper equipment for weighing the aforementioned patients. DETECTO offers innovative designs that can come with mechanical or digital weighing systems. The plastic-molded seats provide needed comfort for patients and all models feature lift-away armrests and footrests.

## 475

### Mechanical Weighbeam Chair Scales

DETECTO's 475 chair scale has a plastic-molded seat for patient comfort complete with lift-away armrests on both sides and a fold-up footrest for easy access during patient transfers from a wheelchair or hospital bed. The die-cast, dual-sided weigh beam is readable from either side and features white enamel finish. No electricity is required, since the scale is mechanical. The 475 has a 450 lb x 4 oz capacity, and the model 4751 is also available with 440 lb x 4 oz / 200 kg x 100 g capacity.



## 6475

### Digital Chair Scales

DETECTO 6475 chair scale provides digital speed and accuracy for weighing while patients are seated in comfort. The 0.7-inch/19-mm high LCD display offer fast readouts accurate to within 0.2 lb / 0.1 kg. The scale's LOCK/RELEASE key holds the patient weight, so that clinical professionals can attend to their patient first before recording the weight. The EMR/EHR-ready RS232 serial port comes standard and Wi-Fi/Bluetooth is available on some models. The scale's capacity is 400 lb x 0.2 lb / 180 kg x 0.1 kg.



## 6880

### Digital Chair Scales with Smart-Phone-Style Indicator

DETECTO's 6880 digital chair scale offers fold-up armrests and footrests to ease patient transfers, 550 lb x 0.2 lb / 250 kg x 0.1 kg capacity, an ergonomically-shaped plastic-molded seat, and smooth gliding 3 in/8 cm diameter rear wheels. Wi-Fi and Bluetooth 5.0 are available on some models for wireless EMR/EHR connectivity. The scale is electronic with a single load cell. Padded hand grips contribute to the 6880 being a pleasure to operate.



### Learn More:

<https://detecto.com/product/product-category/Chair-Scales>

*DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.*

## SOLD BY:

# DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA  
Toll-Free: (800) 641-2008 E-mail: [detecto@cardet.com](mailto:detecto@cardet.com)  
[www.Detecto.com](http://www.Detecto.com)