

# PORTABLE SCALES

DETECTO

MOBILE BARIATRIC WEIGHING

- Spacious weighing platform
- Integral carrying handle
- RS232 serial output
- Rubber non-skid mat
- 1 in / 25 mm high LCD
- Rechargeable built-in battery pack




660lb x 0.5lb 300kg x 0.2kg

Capacity:





## **Spacious weighing platform**






### Easy Set-up and Feet Alignment

The DR660's feet may easily by adjusted for alignment with various weighing surfaces due to transport from one location to another.

## model DR660 PORTABLE SCALES



Capacity	660 lb x 0.5 lb / 300 kg x 0.2 kg
Weight Display	7 digit, 1 in $/$ 25 mm high characters LCD
Dimensions (Display)	10.63 W x 3.54 D x 1.77 H in / 27 W x 9 D x 5 H cm
Dimensions (Platform)	21.5 W x 17.2 D x 1.97 H in / 56 W x 44 D x 5 H cm
Display Cable Length	Flexes from 34 - 126 in / 86 - 320 cm
Power	Built-in rechargeable battery pack (7.2VDC/1200mA) or 120 VAC 60Hz powering a 15 VDC 300 mA wall plug-in UL/CUL listed power supply (included)
Battery Usage	8 hours of continuous use after 14 hours recharging time
Tare	100% of full scale capacity
Temperature	40° to 105°F (5° to 35°C)
Humidity	25% ~95% RH
Keys	HOLD, Ib/kg, ON/ZERO/OFF, TARE
Features	Low battery indicator. Power-saving selectable timing auto shut-off
UPC Code	809161185607
Shipping Weight	28 lb / 12.5 kg

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

### Sold By:

A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA

203 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001 www.Detecto.com