

PH: (417) 673-4631 FAX: (417) 673-5001 WWW.DETECTO.COM

D

ET

ECT

0

## 6855

## High Capacity Waist-High Stand-On Scale

Designed for obese and unsteady patients, this 600 lb capacity scale features a balanced, lightweight design. It's easily tipped back for transportation on its hard, large diameter, rubber wheels yet narrow enough to fit through 30" doors. The battery-powered, LCD display rotates for easy readability. The lowprofile platform, convenient handrail and skid-proof mat add to patient assurance.




Model: Capacity: Platform Size: Display:

**Shipping Weight:** 

**Power:** 

600 lb x 0.2 lb / 270 kg x 0.1 kg 18"W x 14"D / 46 cm x 26 cm • 0.75" high digits, high-contrast LCD • Body Mass Index calculation • LOCK button holds weight after weighing 6 "AA" size batteries (not included) or optional AC adapter. 70 lb / 32 kg



6855

Large LCD and Body Mass Index calculation

75.0

DETECTO