


3300-0001-2A ANALOG BOARD

REVISION RECORD

LTR	DESCRIPTION	DATE	BY
A	RELEASED	2,7,17	28114 bs/
B	CLEANUP	3,24,17	28172 bs/
C	FLOOR SCALE CABLE WAS 3300-0293-1A.	9/19/18	29128 EMJ
D	POWER SUPPLY WAS 6800-1043 POWER CABLE WAS 8555-B428-1A	10,18,18	29164 bs/
E	INDICATOR WAS "ICON", CHNG TO "SOLACE"	1,28,19	29348 bs/



* NOTE: THESE TWO CABLES ARE PART OF 3300-0001-2A ANALOG BOARD ASSEMBLY

UNLESS OTHERWISE SPECIFIED
TOLERANCE ON DIMENSIONS ARE:
 ANGLES ± 1/2°
 INTEGERS/FRACTIONS ± 1/16 IN.
 DECIMALS (.X) ± .03 IN.
 DECIMALS (.XX) ± .01 IN.
 DECIMALS (.XXX) ± .005 IN.
 NOTE: WHOLE NUMBERS MUST BE WRITTEN XX.000 TO INVOKE DECIMAL TOLERANCES.

THIS DRAWING IS THE PROPERTY OF THE CARDINAL SCALE MFG. CO. AND SHALL NOT BE USED NOR REPRODUCED EITHER WHOLLY OR IN PART EXCEPT WITH WRITTEN AUTHORIZATION. ALL RIGHTS OF DESIGN OR INVENTION ARE RESERVED.

Cardinal
 Cardinal Scale Manufacturing Co.
 WEBB CITY, MISSOURI · USA · 64870

CARDINAL DETECTO FULLER WEIGHING SYSTEMS

TITLE DIALYSIS SCALE INTERCONNECT			
SCALE	NTS	DO NOT SCALE DRAWING	NEXT ASSEMBLY
DR.	bs	DATE	B
CH.		DATE	
MODEL		SHEET 1 OF 1	DWG. NO. 3300-0290-02
			REV. E