

EMERGENCY ROOM PRODUCTS

Products you can rely on to help your
ER staff perform at their peak



DETECTO®

AN EMERGENCY ROOM MUST BE READY FOR ANYTHING

DETECTO's medical carts and clinical scales will help you quickly respond to any critical situation with confidence and accuracy.

Medical equipment from DETECTO is built tough to withstand the serious demands of an Emergency Room. Our medical carts are solid, secure, and is customizable with accessories to perfectly fit your needs. You'll be amazed at the variety and flexibility of our scale line, with everything from large roll-on wheelchair scales to touchless physician scales featuring sonar height rod sensors to measure a patient's height and automatically calculate BMI. We also offer waste receptacles and glove box holders in many sizes built to last a lifetime.

RC33669RED-L

**RESCUE Medical Carts
Loaded with Accessories**

DETECTO's steel Rescue medical carts are sturdy and reliable... built to make sure you have the supplies you need where and when you need them. They come in several colors, with red selected for ER use. Carts have your choice of five or six drawers that can be secured with a key or EMG breakaway lock. The ABS countertop is easy to clean and 5-inch wheels feature two parking locks and one steering lock. Optional accessories include IV poles, cardiac boards, catheter holders, drawer dividers, glove dispensers, sharps container holders and waste bins.



8550

**Stretcher Scales with
Column-Mounted Indicator**

DETECTO's Stretcher Scale is perfect for quickly weighing emergency arrivals, with a portable design that makes it easy for anyone to move. The column display is easy to read, with a Hold/Release key to hold the patient's weight after leaving the scale for better patient care. The scale's single-pass weighing feature saves time by allowing you to enter the stretcher or wheelchair tare weight using the numeric keypad on the indicator. You may also use the convenient pushbutton tare feature to remove the stretcher or wheelchair weight if you don't know the weight already. Guide rails and Wi-Fi/Bluetooth available.



6570

**Portable Wheelchair
Scale with Handrail
and Seat**

The 3-in-1 6570 scale provides seated, standing and wheelchair weighing. The fold-away seat and stable handrail offer maximum flexibility without sacrificing portability. The large 32 in x 36 in platform has two-way built-in ramps for maximum accessibility. The scale has a capacity of up to 1,000 lb/450 kg and units are able to be locked in to eliminate transcription errors. Wi-Fi/Bluetooth is available on some models for wireless EMR/EHR connectivity.



APEX-SH

**apex® Digital Physician Scale
with Touchless Sonar
Height Rod**

DETECTO's apex® digital scale with sonar height rod lets you measure height and weight without pushing a single button. The large .75-inch blue LCD readout is easy to view. The scale provides BMI calculation, runs on AC or battery power, and comes with a wide variety of connectivity options including RS232 serial and micro USB, with some models featuring Wi-Fi/Bluetooth and Welch Allyn connectivity for electronic medical records.



P-100R

**Step-On Cans, the Rugged
Waste Receptacle**

Named for their rugged toughness, Step-On Cans are perfect for hospitals, doctors' offices, patients' rooms, and laboratory settings. Equipped with a heavy-duty foot pedal and linkage and a bag securing mechanism, they come in a variety of five different sizes and capacities (4, 6, 8, 12 and 25 gallon) and two different materials: stainless steel or a baked epoxy finish (available in red or white). The cans are flame retardant, have a self-closing lid for sanitation, and feature unitized construction.



GH3

**Universal
Wall-Mount Steel
Glove Box Holders**

Conserve counter space with DETECTO's universal glove box holders. The open design allows you to see what size or type of glove is being dispensed. The holders accommodate most glove boxes and are universal enough to be used for tissue holders where counter or shelf space is limited. Models capable of holding one, two, or three boxes of gloves. Lifetime warranty.



DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

SOLD BY:

DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com