


LONG-TERM CARE PRODUCTS

DETECTO manufactures the safest and most accessible line of weighing products for long-term care patients when it comes to wheelchair scales, chair scales, and in-bed scales.



DETECTO®

LONG-TERM CARE DESERVES LONG-LASTING DURABILITY

DETECTO products are mainstays in long-term care facilities and nursing home residences for a number of reasons. DETECTO brand products are known for quality and safety with long-lasting dependable accuracy. The DETECTO line includes wheelchair, chair, in-bed, and eye-level scales designed to comfortably accommodate frail, elderly patients who need assistance while weighing. DETECTO also offers plastic and metal waste receptacles and medical carts for waste disposal and mobile storage.

6550

Column-Mounted Weight Display with Collapsible Column for Transport



DETECTO's 6550 portable wheelchair scales offer time-saving and patient comfort benefits. The 6550 features 1,000 lb / 450 kg capacity, Wi-Fi/Bluetooth on some models, vertical fold-up storage, two-way built-in ramps for accessibility from both sides, and 32 in x 36 in / 81 cm x 91 cm platform. With only one hand you can quickly slide the release lever and fold down the 6550's column down for easy transport or vertical storage.

7550

Wall-Mount Fold-Up Scale



The 7550 wall-mount wheelchair scale features a fold-up design that saves space when not in use. It can be installed in high-traffic areas for convenient weighing, and fold up when finished weighing. The scale accommodates most wheelchairs, and yet requires very little floor space when folded up. The 7550 features a full-color, touch-screen display and Wi-Fi/Bluetooth for wireless connectivity with patient record systems.

6570

Portable Wheelchair Scale with Handrail and Seat



The 3-in-1 6570 scale provides seated, standing and wheelchair weighing. The fold-away seat and stable handrail offer maximum flexibility without sacrificing portability. The large 32 in x 36 in platform has two-way built-in ramps for maximum accessibility and 1,000 lb/450 kg capacity.

BRW1000

Larger-Size 32 in x 40 in Platform for Extra- Wide Bariatric Wheelchairs



DETECTO's BRW1000 is designed to accommodate extra-wide bariatric wheelchairs with its large 32 in x 40 in x 2.2 in / 81 cm x 101 cm x 5.6 cm platform. Since the display is connected to the scale by a long coiled cable that extends up to 14 ft / 4.2 m, it can be positioned where the healthcare professional can most easily read the results. Models are available with Wi-Fi/Bluetooth and also Welch Allyn connectivity. The BRW1000 can accommodate motorized wheelchairs.

6475

Digital Chair Scale with 400 lb Capacity



The ultimate in patient safety and comfort, DETECTO's digital chair scales are completely mobile, with built-in wheels and battery power. Lift-away armrests and a fold-up footrest make it easy for transferring dialysis patients, disabled patients, and the elderly. RS232 included standard and Wi-Fi/Bluetooth on some models.

6880

Chair Scale with Liftaway Armrests



DETECTO's 6880 chair scale is the simple solution for patients with decreased mobility and features 550 lb / 250 kg capacity, dual wheels for transport, smartphone-style indicator, liftaway padded armrests and footrests, and battery or AC power. Wi-Fi and Bluetooth wireless connectivity models are available.

RC33669WHT

RESCUE Steel Medical Carts with Customizable Accessories



DETECTO's Rescue medical carts are sturdy and reliable... built to make sure you have the supplies you need where and when you need them. Carts have your choice of five or six drawers that can be secured with a key or EMG breakaway lock and the ABS countertop is easy to clean.

PL600

Patient Lift Scale



DETECTO's Patient Lift Scale allows for a person to be suspended and weighed from a portable lift mechanism. The lock-in weight feature stores the last reading whenever the unit is turned off. PL600 provides easy setup and operational use.

P-100

Step-On Cans, the Rugged Waste Receptacle



Step-On Cans come equipped with a heavy-duty foot pedal and linkage and a bag securing mechanism, they come in a variety of five different sizes and capacities and two different materials: stainless steel or a baked epoxy finish (available in red or white).

WM-18RED

Waste Mate Plastic Waste Receptacles



Available in four sizes and red and beige colors, DETECTO's Waste Mate commercial-grade plastic waste receptacles feature an easy-to-use foot pedal that flips open the lid for no-touch trash disposal. The integrated internal clamps secure the trash bag to the rim of the can to prevent it from slipping into the can.

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

SOLD BY:

DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com