

Baby Scale



6745 Digital Baby Scale

- 30 lb / 15 kg capacity with 1 gram divisions
- 23" x 9" platform (flat portion)
- Touchscreen LCD and Wi-Fi/Bluetooth included

Physician Scales



APEX Series Column Scales with Mechanical Height Rod

- 600 lb x 0.1 lb / 270 kg x 0.05 kg
- 17" x 17" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



APEX-RI Series Low-Profile Physician Scales

- 600 lb x 0.1 lb / 270 kg x 0.05 kg
- 17" x 17" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



APEX-SH Series Column Scales with Touchless Sonar Height Rod

- 600 lb x 0.1 lb / 270 kg x 0.05 kg
- 17" x 17" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



ICON Series Column Scales with Sonar Height Rod

- 800 lb / 400 kg capacity with 0.1 lb / 0.05 kg divisions
- 17" x 17" platform
- Touchscreen LCD and Wi-Fi/Bluetooth included

Chair Scale



6880 Series Rolling Chair Scales

- 500 lb x 0.2 lb / 250 kg x 0.1 kg
- Foldaway armrests and footrest
- Available with AC adapter and/or Wi-Fi/Bluetooth

Wheelchair Scales



6400 Series Wheelchair Scales with Basic-Level Indicator

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 32" x 36" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth




6500 Series Wheelchair Scales with Upgraded Indicator

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 32" x 36" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6600 Series Bariatric Wheelchair Scales with Basic-Level Indicator

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 32" x 40" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



BRW1000 Series Bariatric Wheelchair Scales with Upgraded Indicator

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 32" x 40" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6550 Series Wheelchair Scales with Column-Mounted Indicator

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 32" x 36" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6560 Series Wheelchair or Standing Scale with Curved Handrail

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 32" x 36" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6570 Series Wheelchair, Seated, or Standing Scale with Handrail

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg
- 32" x 36" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



NTEP Certified Medical Scales

The Only Medical Scales with Verified Accuracy

Bariatric Scales



6800 Series

Low-Profile Bariatric Scale

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 24" x 24" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6855 Series

Handrail Scale

- 600 lb x 0.2 lb / 270 kg x 0.1 kg capacity
- 18" x 14" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6855MHR Series

Handrail Scale with Mechanical Height Rod

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 18" x 14" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6854DHR Series

Handrail Scale with Digital Height Rod

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 18" x 14" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6856 Series

Handrail Scale

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 24" x 24" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6856MHR Series

Handrail Scale with Mechanical Height Rod

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 24" x 24" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6857DHR Series

Handrail Scale with Digital Height Rod

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 24" x 24" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth



6868 Series

Bariatric Flip Seat Scale

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 24" x 24" platform
- Available with AC adapter and/or Wi-Fi/Bluetooth

Wet Diaper Scale



AP-4KD

Compact Digital Countertop Scale

- 4,000 g x 1 g
- 9" x 6" stainless steel tray
- Comes with AC adapter
- 0.5" red LED display



DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA
P: 417-673-4631 or 1-800-641-2008 F: 417-673-2153
www.DETECTO.com