



# PEDIATRIC CART

## 9 Multi-Color Drawers



*DETECTO model  
RC33333336PED*

DETECTO's Rescue series Pediatric cart is designed for use with the proven clinical methodology for treating life-threatening emergencies in pediatric patients. When seconds count, the color-coded drawers enable doctors and clinicians to quickly identify critical supplies while responding to pediatric emergencies.

Contents are secured with individual plastic seals for each specific drawer or a keyed lock which provides immediate access to all drawers and emergency supplies.



# DETECTO Rescue Series Pediatric Medical Cart with Push Handle and User-Friendly Design

DETECTO's Rescue series pediatric carts come in 9 multi-colored drawers using the Broselow® color-coding method and have 5-inch diameter wheels which feature two parking locks and one steering lock. Rescue carts are manufactured from quality steel cart body construction that is powder-coat painted for durability with ABS countertops that are easy to clean. ABS bumpers protect all four corners of the medical carts. The push handle offers a user-friendly design.



- White frame with 9 multicolor drawers
- Drawer Sizes: 3 in cream / 3 in pink/red / 3 in lavender / 3 in yellow / 3 in white / 3 in light blue / 3 in bright orange / 3 in fresh green / 6 in cream
- Each drawer has an individual breakaway lock
- Keyed lock for all drawers
- Numerous optional accessories available
- Push handle with user-friendly design
- Auto-return rail on drawers
- ABS countertop surface is easy to clean
- Quality steel cart body construction
- 2 parking locks and 1 steering lock on casters
- 5-in diameter wheels
- ABS bumpers on each of the 4 corners



Individual drawer locks and universal keyed lock for all drawers. Optional CAPS red plastic seals are available that come in qty. 50, individually numbered.

SPECIFICATIONS	
<b>Model</b>	RC33333336PED
<b>Overall Dimensions (without push handle)</b>	26.9 in W x 21.5 in D x 45 in H / 68.3 cm W x 54.6 cm D x 114.3 cm H
<b>Overall Dimensions (with push handle)</b>	29.8 in W x 21.5 in D x 47 in H / 75.7 cm W x 54.6 cm D x 119.4 cm H
<b>Inner Drawer Dimensions (3-inch-high drawers)</b>	22.3 in W x 15.5 in D x 2.5 in H / 56.6 cm W x 39.4 cm D x 6.4 cm H
<b>Inner Drawer Dimensions (6-inch-high drawers)</b>	22.3 in W x 15.5 in D x 5.5 in H / 56.6 cm W x 39.4 cm D x 14.0 cm H
<b>Shipping Weight</b>	161 lb / 73 kg
<b>Shipping Dimensions</b>	28 in x 22 in x 46 in / 71 cm x 56 cm x 117 cm
<b>Options</b>	CAPS red plastic seals (qty. 50, individually wrapped) and a wide range of optional accessories including bin organizers, cardiac board, defibrillator shelf, IV pole, drawer dividers, sharps container holder, catheter holder, oxygen tank holder, and waste bin.

**SOLD BY:**

 **DETECTO**  
www.Detecto.com

203 E. Daugherty, Webb City, MO 64870 USA  
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153