

# PEDIATRIC

scale cart



# PSCART



Make every pediatric weigh-in easier with the DETECTO Pediatric Scale Cart. Purpose-built for the 6745 scale, this mobile cart provides the perfect working height for medical staff while keeping essential supplies close at hand.

Designed with safety and efficiency in mind, the cart's protective lip prevents items from sliding off, while smooth-rolling casters allow quiet, effortless movement throughout patient rooms.

With added storage for diapers, wipes, and other essentials, the Pediatric Scale Cart keeps everything you need within reach - so you can focus on patient care.



# DETECTO®



- **Optimized Height:** Provides the ideal working height for nurses when weighing and measuring infants, reducing strain and improving ergonomics.
- **Integrated Storage:** Enclosed compartment with fixed shelf holds diapers, wipes, bandages, and other frequently used medical supplies for greater efficiency.
- **Custom Fit for the DETECTO 6745 Scale:** Specifically designed to ensure stability and reliability during use.
- **Protective Edge Lip:** Keeps the scale and other items securely in place, preventing accidental sliding or movement.
- **Smooth Mobility:** Quiet, hospital-grade casters glide easily across rooms and hallways, with secure stability when locked.
- **Compact & Efficient:** A space-saving solution that combines scale positioning and supply storage in one mobile cart.

SPECIFICATIONS	
<b>Model</b>	PSCART
<b>Construction</b>	Quality aluminum cart body with powder-coat paint and ABS countertop and bumpers
<b>Wheels</b>	4-inch-diameter wheels (2 parking locks and 1 directional steering lock)
<b>Cart Dimensions (w x d x h)</b>	33.07 in x 21.63 in x 30.94 in / 84.01 cm x 54.94 cm x 78.58 cm
<b>Cabinet Dimensions (w x d x h)</b>	28 in x 18.5 in x 21.5 in / 71.12 cm x 46.99 cm x 54.61 cm
<b>Shelf Height</b>	Bottom – 9.875 in / 25.08 cm, Top – 10.875 in / 27.62 cm
<b>Shipping Weight</b>	77 lb
<b>Maximum Loaded Weight of Cart</b>	1,100 lb / 498 kg
<b>Warranty</b>	12-year limited warranty

