

PHYSICIAN SCALES

From state-of-the-art sonar height rods and advanced electronics to traditional weighbeam physician scales



DETECTO

EYE-LEVEL WEIGHT DISPLAY CONVENIENCE OR LOW-PROFILE PORTABILITY

DETECTO is the largest medical scale manufacturer in the US, and the flagship product from the beginning has been the company's well-known physician scale. DETECTO's traditional weigh beam eye-level design has been around for over 120 years, being used by doctors, nurses, and healthcare professionals around the world. DETECTO has also developed state-of-the-art digital physician scale models to meet the demands of today's medical world, including touchless sonar height rods.

solo®

Eye-Level Digital Scale with Inline Mechanical Height Rod

DETECTO's solo® eye-level clinical scale is the most economical, high-quality digital physician scale with height rod available. The solo® scale features 550 lb x 0.2 lb / 250 kg x 0.1 kg capacity, 14 in W x 15 in D x 2.4 in H low-profile platform, inline mechanical height rod (measures 30 in to 84.5 in / 76 cm to 214.5 cm), and a RS232 serial port.



apex® Series

Eye-Level Digital Scales with Inline Mechanical or Sonar Height Rods

DETECTO's apex® eye-level digital physician scales feature inline mechanical height rods or touchless sonar height rods and an extra-wide, flat platform measuring 17 in W x 17 in D x 2.75 in H. The apex® scale features 600 lb x 0.2 lb / 300 kg x 0.1 kg capacity and 0.75-inch-high clinical-blue LCD weight readouts. Wi-Fi/Bluetooth models are available.



icon®

Premium Eye-Level Digital Scale with Touchless Sonar Height Rod

DETECTO's icon® eye-level digital physician scales feature state-of-the-art sonar height rods and an ultra-thin platform only 1.5 inches high. The icon® series scales feature an extra-high 1,000-lb capacity, multi-color touchscreen LCD display, 0.9-inch-high LCD weight digits, Wi-Fi/Bluetooth for wireless EMR/EHR, and a 17 in W x 17 in D weighing platform.



Weighbeam Series

Eye-Level Mechanical Weighbeam Scales With Height Rod, Handpost, and/or Wheels


DETECTO's world-renowned weighbeam physician scales feature a 10.5 in x 14.5 in (27 cm x 37 cm) platform and dual-reading die-cast weighbeam which may be read from either side of the scale. Thirteen different models are available in lb, lb/kg, and kg with capacities of 440 lb/200 kg, 450 lb, or 200 kg and accessories such as handposts, height rods, and wheels.



APEX-RI Series

Digital Low-Profile Scales with Remote Display

DETECTO's APEX-RI series portable scales with remote indicators feature an extra-wide platform measuring 17 in W x 17 in D x 2.75 in H, 600 lb x 0.2 lb / 300 kg x 0.1 kg capacity, 0.75-inch-high clinical-blue LCD weight readouts, 6 ft/1.8 m cable from the scale base to indicator, and integral carrying handle. Wi-Fi/Bluetooth models are available.



SOLO-RI

Digital Low-Profile Scale with Remote Display

DETECTO's SOLO-RI scale has a low-profile scale base with 7.25 ft / 2.2 m cable to the remote display, so you can position the indicator wherever it's most convenient for viewing. The scale features 550 lb x 0.2 lb / 250 kg x 0.1 kg capacity, 14 in W x 15 in D x 2.5 in H platform, Body Mass Index, 1-inch-high LCD digits, 6 AA batteries (included), OneWeigh units locking, and countertop or wall-mount bracket (included).



Learn More:

<https://detecto.com/product/product-category/Physician-Scales>

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

SOLD BY:

DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com