



# ProDoc® PD SERIES PROFESSIONAL DOCTOR SCALES

MEDICAL GRADE ACCURACY: WEIGHT • HEIGHT • BODY MASS INDEX



DETECTO  
Model PD100



- 480 lb x 0.2 lb / 220 kg x 0.1 kg
- Space saving, compact design
- 4-user memories and recall function
- 6 AA batteries (not included) or optional AC adapter
- Solid, stable base for patient assurance and accuracy
- Body Mass Index (BMI) calculation
- Metric conversion (Weight and Height)
- Selectable Sleep Mode (up to 9 minutes)
- Selectable Auto Shut-off Feature (up to 9 minutes)
- Large 1-inch / 25-mm high LCD readouts



DETECTO  
Model PD200



# PD SERIES PROFESSIONAL DOCTOR SCALES



## Model PD100 Low-Profile Floor Scale with Remote Display

The PD100's remote display may be positioned on desks, counters, or even mounted on the wall (shown below) up to 87 inches / 220 cm from the scale. The LCD provides easy-to-view 1-inch-high (25 mm) readouts.

Optional model PRODOC-CASE carrying case for PD100



Ready for weighing use right out of the box (No assembly required).



## Model PD200 Space-Saving Physician Scale

The PD200 features a low-profile, angled display for easy measurement viewing. The PD100 and PD200 provide BMI (Body Mass Index) calculation by conveniently selecting your height with easy-to-use scroll keys while weighing.



### PD SERIES SCALE SPECIFICATIONS

<b>Model Number</b>	PD100 (Floor Scale with Remote Display) PD200 (Floor Scale)
<b>Capacity</b>	480 lb x 0.2 lb / 220 kg x 0.1 kg
<b>Weight Units</b>	Pounds or Kilograms (selectable)
<b>Power Requirements</b>	Six (6) "AA" size alkaline batteries (not included) OR an optional Medical device 9V AC/DC wall plug-in adapter (Cardinal part number PD-AC, includes USA plug). Also available UK plug (PD-UKPLUG) or EU plug (PD-EUPLUG).
<b>Display</b>	Five digit, seven segment, 1.0 in (25 mm) high LCD
<b>Operating Temperature</b>	50 to 104° F (+10 to +40° C)
<b>Function Keys</b>	ON/OFF, BMI HEIGHT/ENTER, LOCK/RELEASE, UNITS, ZERO, RECALL and STORE
<b>Dimensions</b>	
PD100 Display	5.25 in H x 7.87 in W x 2.125 in D (134 mm x 200 mm x 54 mm)
PD100 Base	3.33 in H x 12.5 in W x 13.0 in D ( 85 mm x 318 mm x 330 mm)
PD200	8.75 in H x 12.6 in W x 16.9 in D (222 mm x 320 mm x 430 mm)

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.



203 E. Daugherty, Webb City, MO 64870 USA  
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001  
www.Detecto.com

Sold By:

