



**Slim PRO**

**LOW-PROFILE HEALTHCARE SCALES**

# Slim PRO

440 lb x 0.2 lb / 200 kg x 100 g



\* shown with optional serial interface cable and model P50 direct thermal printer

## LOW-PROFILE HEALTHCARE SCALES

*...a step ahead of the competition*

### BATTERY POWER

The SlimPRO is powered by 4 AAA batteries, which are included with the scale. Auto-On and Auto-Off functions help conserve battery life.

### LARGE LCD READOUTS

The SlimPRO's big 1.5 in / 3.8 cm high LCD digits are easily viewable. Convenient units switching is performed through a switch on the back of the scale.

### SWIVEL FEET

DETECTO's SlimPRO can even weigh accurately on carpets. The special swivel feet allow the scale to be used on a variety of floor surfaces.

- Integral carrying handle for convenience.
- Incredible 1.4-in/3-cm thin construction makes stepping on and off very easy.
- Extra-wide 17-in/32-cm weighing platform accommodates patients up to 440 lb/200 kg.
- RS232 serial port allows connectivity with PC for EMR/EHR or optional Detecto P50 printer.
- Designed for performance and durability in clinical-grade usage weighing heavy patients.
- The scale automatically turns on when a patient steps on the scale and turns off when they step off.
- Weight is conveniently held onscreen for 5 seconds after weighing.

