



Slim PRO

TALKING
LOW-PROFILE
HEALTHCARE
SCALES



- 600 lb / 272 kg Capacity
- Backlit Blue LCD
- Bluetooth Connectivity
- Mother/Baby Weighing Mode
- Talking Scale: Speaks 6 Different Languages

CLEAR PLEASANT VOICE ANNOUNCEMENT OF WEIGHT

"You weigh 127.9 pounds."

Easy to use for the vision impaired

DETECTO's SlimPRO talking healthcare scales are ideal for the elderly, obese, blind, and telemedicine use. The scale is easy to use for the vision impaired and provides a clear, pleasant voice announcement of the weight in your choice of six languages: English, Spanish, French, Arabic, Hindu, or Chinese (switch selectable). The voice function can also be turned off.

The SlimPRO is a high-capacity, low-profile digital scale for clinical or personal weighing use. The spacious, extra-wide weighing platform easily accommodates patients up to 600 lb/272 kg. The stylish gray color blends into any environment. The scale automatically turns on when a patient steps on the scale and turns off when they step off.



THIN, LOW-PROFILE HEALTHCARE WEIGHING

As the name implies, DETECTO's SlimPRO is an incredibly thin 1.4-in/3-cm high scale that contains loads of weighing performance and durability for clinical-grade usage weighing heavy patients.

- Extra-wide platform for patients of all sizes
- Large 1.5 in/38 mm high LCD digits for easy reading
- Auto off function conserves battery life
- Body Mass Index calculation
- 4 AAA batteries included
- Step-On/Off function for weighing without switching on scale
- Ribbed black mat for weighing traction and user safety



BLUETOOTH CONNECTIVITY

The SlimPRO features Bluetooth Low Energy (BLE) 4.0 for wireless connectivity with mobile apps.



MOTHER/ BABY MODE

Mother/Baby mode factors the weight of infants and toddlers held by an adult.



LARGE BACKLIT BLUE LCD FOR EASY VISIBILITY

Thin profile for easy accessibility combined with high 600-lb/ 272-kg capacity

- Non-skid rubber mat provides safety when stepping on and off the scale
- Low battery indicator for added convenience
- Flat construction makes stepping on and off very easy
- Weight is conveniently held onscreen for 5 seconds after weighing
- Integral carrying handle built into the scale
- Sturdy, well-built construction for long-term durability
- No need to turn it On and Off - the scale does it for you

Detecto Model SlimPRO
Digital Low-Profile
Floor Scale



IDEAL FOR:

- Elderly
- Obese
- Blind
- Telemedicine



DETECTO
www.Detecto.com



LARGE LCD READOUTS

The SlimPRO's big 1.4 in / 36 mm high LCD digits are easily viewable. Convenient units switching is performed through a switch on the back of the scale.

SWIVEL FEET

DETECTO's SlimPRO can even weigh accurately on carpets. The special swivel feet allow the scale to be used on a variety of floor surfaces.

BATTERY POWER

The SlimPRO is powered by 4 AAA batteries, which are included with the scale. Auto-On and Auto-Off functions help conserve battery life.

The SlimPRO may be used in applications including physician offices, clinics, gyms, home healthcare, and bathroom personal use.





Slim PRO



TALKING LOW-PROFILE HEALTHCARE SCALES

Model	SLIMPRO	
Capacity	600 lb x 0.1 lb / 272 kg x 0.1 kg	
Display Type	Four-digit, seven segment LCD, 1.4 in / 36 mm high digits with blue backlight	
Display Area Size	3.15 in W x 1.6 in H / 8 cm W x 4 cm H	
Platform Dimensions	17 in W x 14.5 in D x 1.4 in H / 43 cm W x 37 cm D x 3 cm H	
Platform Surface	Ribbed black mat for weighing traction and user safety	
Voice	Talking Scale Speaks 6 Languages: English, Spanish, French, Arabic, Hindu, Chinese, OFF	
Indicators	Weight, Height, BMI Calculation, Low Battery, Mother and Baby Mode	
Connectivity	Bluetooth Low Energy (BLE) 4.0	
Power	Four (4) "AAA" Batteries (included)	
Enclosure Construction	Polycarbonate	
Units	Kilograms, Pounds	
Functions	Mother and Baby weighing mode	Baby weight starts from 2.2 lb (1 kg)
	Auto-ON	Requires minimum weight of 11 lb (5 kg)
	Auto-OFF	Approximately 5 seconds
	Unit	lb/kg, selected by switch on bottom of scale
	BMI Calculation	Manual Input Height (75 cm to 225 cm with 1 cm)
Included	4 carpet feet, 4 AAA batteries, and product manual	
Shipping Weight	5 lb / 2.3 kg	
UPC Code	809161189001	

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

 **DETECTO**
www.Detecto.com

102 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153

SOLD BY: