

DETECTO®

Slim TALK



Your weight is one-hundred eighty point six pounds

TALKING SCALES (BILINGUAL ENGLISH/SPANISH)



Slim TALK XL



Su peso es setenta y ocho punto nueve kilogramos

IDEAL FOR:

- Elderly
- Obese
- Blind
- Telemedicine

TALKING SCALES (ENGLISH AND SPANISH)

Easy to use for the vision-impaired

- Two capacities:
400 lb x 0.1 lb / 180 kg x 0.1 kg (SlimTalk)
or 550 lb x 0.1 lb / 250 kg x 0.1 kg (SlimTalk XL)
- Extra-wide platform for patients of all sizes
- 5 carpet feet included
- Large 1.5 in/38 mm high LCD digits for easy reading
- 3 AAA batteries included
- Auto off function conserves battery life
- Sturdy construction for long-term durability
- Clear, pleasant voice announcement of weight
- Very low-profile platform only 1.4 in/36 mm high for easy access
- Step-On/Off function for weighing without switching on scales



NON-SLIP TEXTURED SURFACES



Plastic platform with textured surface for weighing traction and user safety.



Plastic platform with black rubber mat for weighing traction and user safety.

Model	SLIMTALK	SLIMTALKXL
Capacity	400 lb x 0.1 lb / 180 kg x 0.1 kg	550 lb x 0.1 lb / 250 kg x 0.1 kg
Languages	Two (English and Spanish)	
Overall Dimensions	14.3 in W x 12.3 in D x 1.4 in H / 363 mm W x 312 mm D x 36 mm H	
Display Size	3.2 in W x 1.7 in H / 81 mm W x 43 mm H	
Digit Height	1.5 in / 38 mm	
Weight Display	Four digit, seven segment LCD with 1.5 in / 38 mm high digits	
Voice	English / Spanish / Off	
Power	3 "AAA" batteries included (4.5 VDC)	
Indicators	lb, kg, Low Battery	
Functions	Auto-ON (minimum weight of 11 lb / 5 kg is required) Auto-OFF (10 seconds)	
Platform Construction	Plastic	
Included	5 carpet feet, 3 AAA batteries, and product manual	
Net Weight	5.7 lb / 2.6 kg	
Shipping Weight	7 lb / 3.2 kg	
UPC Code	809161302608	809161302707

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

SOLD BY:



102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com