

# DETECTO®

## solace

IN-FLOOR DIALYSIS SCALES



## IN-FLOOR DIALYSIS PATIENT WEIGHING

- Flush-Mounted with Floor for Easy Accessibility
- Stainless Steel, Anti-Corrosion Weighing Platform
- Multi-Color, Touchscreen Indicator/Printer Combo Included



DETECTO's Solace series are available in 6 x 4 ft sizes for weighing stretchers.



DETECTO's Solace series flush-mounted platform scales eliminate tripping hazards and hallway accessibility issues.



**solace**  
IN-FLOOR DIALYSIS SCALES



# IN-FLOOR DIALYSIS SCALES

DETECTO's Solace series in-floor dialysis scales are flush-mounted platform scales that feature stainless steel decks and include a pit frame and trim ring. The flush surface eliminates tripping hazards and hallway accessibility issues.

The Solace series is available with or without a handrail and all models include a recessed-mount indicator and printer combination. These in-floor scales easily accommodate patients with walkers, canes, wheelchairs, or seated in dialysis chairs.



*The ID models are designed to be mounted in a pit flush with the floor. A tubular handrail with padded handgrip is available on some models for walk-on patients.*

## DIALYSIS + BARIATRIC WEIGHING

- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- 8 standard models to choose from
- Standard platform sizes: 3 x 3 ft, 4 x 3 ft, 4 x 4 ft, and 6 x 4 ft
- Handrails available for patient stability while weighing
- Indicator/printer, pit frame, scale, and trim ring included
- Custom sizes and configurations are available



# RECESSED-MOUNT INDICATOR + PRINTER COMBO INCLUDED



The integrated scale displays on DETECTO's multi-color indicator that is wall-mounted with a stainless steel cover making it easy for clinicians to read and accurately record measurements.

DETECTO's Solace series in-floor scales are covered with durable stainless steel to provide years of maintenance-free use. The built-in thermal printer makes it easy to print out a weight record for patient recordkeeping.

- Multi-color, touchscreen LCD display
- Built-in thermal ticket printer (included)
- Displays weight, height, and BMI onscreen simultaneously
- Large 0.9-inch-high, clinical-blue LCD weight digits
- Push button or manual entry of tare weights for wheelchairs and stretchers
- Automatically locks on weight temporarily for a few seconds in order to easily view and record measurement
- Menu button provides quick access to user-friendly menu settings for scale customization
- Standard 110V AC power 110V – 240V AC universal switching power supply

ID: 5432  
GROSS: 489.0 lb  
TARE: 100.0 lb  
NET: 389.0 lb

ID: 8976  
GROSS: 218.2 lb  
TARE: 50.0 lb  
NET: 168.2 lb

*The Solace series in-floor scale's ticket printer can print Gross, Tare, Net and Patient IDs for record keeping.*



The Solace series is flush mounted with the floor and includes the pit frame and trim ring, plus an interface cable from the scale to indicator.



DETECTO's Solace series in-floor scales are ideal for dialysis clinics, emergency rooms, and clinical facilities where fast, efficient patient measurement is necessary.

## Complete in-floor scale package for continuous patient weight monitoring



In-floor scales comes with model 855RMP digital weight indicator/printer combo, pit frame, scale, stainless steel deck, and trim ring. Custom sizes and configurations are available.

## SOLACE SERIES IN-FLOOR DIALYSIS SCALES



Model	ID-3636S-855RMP	ID-3636SH-855RMP	ID-4836S-855RMP	ID-4836SH-855RMP
<b>Capacity</b>	1,000 lb x .2 lb / 450 kg x .1 kg			
<b>Platform Size</b>	36 x 36 in / 91 x 91 cm	36 x 36 in / 91 x 91 cm	48 x 36 in / 122 x 91 cm	48 x 36 in / 122 x 91 cm
<b>Handrail</b>		●		●
<b>Power</b>	110V – 240V AC universal switching power supply (included)			
<b>Shipping Weight</b>	255 lb / 115 kg	265 lb / 120 kg	311 lb / 141 kg	321 lb / 146 kg
<b>UPC Code</b>	809161203707	809161203806	809161203905	809161204001

## SOLACE SERIES IN-FLOOR DIALYSIS SCALES



Model	ID-4848S-855RMP	ID-4848SH-855RMP	ID-7248S-855RMP	ID-7248SH-855RMP
<b>Capacity</b>	1,000 lb x .2 lb / 450 kg x .1 kg			
<b>Platform Size</b>	48 x 48 in / 122 x 122 cm	48 x 48 in / 122 x 122 cm	72 x 48 in / 183 x 122 cm	72 x 48 in / 183 x 122 cm
<b>Handrail</b>		●		●
<b>Power</b>	110V – 240V AC universal switching power supply (included)			
<b>Shipping Weight</b>	388 lb / 176 kg	398 lb / 181 kg	512 lb / 232 kg	522 lb / 237 kg
<b>UPC Code</b>	809161204100	809161204209	809161204308	809161204407

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.



102 East Daugherty St., Webb City, MO 64870 USA  
 Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com  
 www.Detecto.com

**SOLD BY:**