

SUPPORTING THE TINIEST LIVES IN THE NICU

At DETECTO, we understand that having a baby in the NICU is a stressful and scary experience. That is why we are committed to providing precise, reliable scales to help alleviate that stress and empower caregivers with confidence.



DIGITAL BABY SCALE



6745

- 0 – 10 lb x 0.002 lb / 10 – 30 lb x 0.005 lb
- 0 – 10 lb x 0.02 oz / 10 – 30 lb x 0.05 oz
- 0 – 5 kg x 0.001 kg / 5 – 15 kg x 0.002 kg
- 0 – 5,000 g x 1 g / 5,000 g – 15,000 g x 2 g
- Full-color, touchscreen 3.2 in (diagonal) display
- Roomy 22 in x 9 in / 57 cm x 23 cm weighing tray
- Integral measuring tape (0-22 in / 0-56 cm)
- Wi-Fi and Bluetooth 5.0 for EMR/EHR connectivity

WET DIAPER SCALES



WPS12DT

- 12 lb x 0.1 oz / 12 lb x 0.005 lb / 12 lb x 1/8 oz / 192 oz x 0.05 oz / 192 oz x 1/8 oz / 5,500 g x 1 g
- Fully submersible underwater for easy washdown
- Stainless steel watertight scale and lift-off tray



AP-4KD

- 4,000 g x 1 g
- 9 in x 6 in x 2 in weighing tray
- Stainless steel scale and lift-off tray
- Easy-to-read, 0.5-in-high red digital LED display

DETECTO®