## DETECTO VET-70 DIGITAL VETERINARY SCALE

Precision, Portability, Performance for Weighing Small Animals



DETECTO's model VET-70 compact scale features a stainless steel platform for long-lasting durability and easy wipe-down cleaning after animal weighing. The scale weighs in .02 lb / .01 kg increments, yielding precision results for small animals, such as cats, dogs, and rabbits. The plug-in AC adapter and lightweight build allow for easy set-up and mobility. DETECTO's unique StableSense digital filtering greatly reduces weight fluctuations resulting from animal movement, allowing for correct weights when dealing with active pets.

- Platform Size: 12 x 14 in / 31 x 36 cm
- Power: 100-240 VAC 50/60Hz 12 VDC 1A wall plug-in UL/CSA listed AC power adapter included (Cardinal part no. 6800-1045)
- Display: 5-digit, 7-segment, 0.875-in / 22-mm high backlit bright green LCD with 0.71-in / 18-mm high digits
- Keys: ZERO/POWER, UNITS
- Construction: Stainless steel weighing platform, ABS plastic indicator and mounting bracket, and steel base
- Shipping Weight: 26 lb / 11.8 kg





Ph. (800) 641-2008 • detecto@cardet.com 203 East Daugherty St., Webb City, MO 64870