

- 330 lb x 0.2 lb / 150 kg x 100 g Capacity
- 36 in W x 22 in D / 90 cm W x 55 cm D Stainless Steel Platform
- Rechargeable Battery Power and AC Adapter Included

# ALL-PURPOSE DIGITAL SCALE

## Comfortable Weighing FOR All Size Animals

DETECTO's economical, portable VET-330WH digital scale with remote weight display is perfect for weighing small to large size animals. DETECTO's VET-330WH digital platform scale features rechargeable battery power and comes with an AC adapter. The scale's easy-to-read remote LCD weight display features 1-inch-high digits that are bold enough to be easily read across the room. The VET-330WH features a low battery indicator and power-saving selectable timing auto shut-off. The stainless steel platform provides long-lasting durability and easy mobility and versatility.

The VET330WH's rechargeable battery power combined with its light weight and integral wheels/carrying handle allow the scale to be used in numerous locations. The HOLD key on the remote display allows the animal's weight to be retained after leaving the platform. The remote display's lengthy 7 ft/2.1 m coiled cable allows it to be positioned away from the platform or nearby, wherever weight display is convenient for viewing.



#### Featuring Integral Wheels, Guide Handle, and Rechargeable Batteries for Portable Weighing

- TARE key (combined with 330 lb capacity) allows you to hold the animal in your arms while weighing.
- Mount the display permanently on the wall or temporarily on any countertop using the universal bracket.
- Includes four leveling feet to accommodate slightly uneven floor surfaces.

° 2 190. | 0000

- Stainless steel, scratch-resistant platform for durability and easy wipe-down cleaning.
- 7.0 ft (2.1 m) cable from scale base to display allows the indicator to be placed wherever it's easiest to operate and view.
- Included wall or desk mount display bracket and mounting hardware.
- Minimal assembly required (simply install non-skid feet and mounting bracket to display).
- Rechargeable battery included with low battery indicator.
- Includes an AC power adapter to power the scale and charge its battery.
- Selectable Auto Shut-Off feature (up to 5 minutes) turns the scale off when not used which extends battery life.
- Zero is established on power up routine and maintained by auto-zero circuitry.
- HOLD key allows the animal's weight to be retained after leaving the platform.





The scale's large-size platform allows it to be used for a wide range of weighing purposes: animal health, commercial laundry, feed and grain, food distribution, meat processing, and brewery distribution.



The scale comes with a UL/CUL listed AC power adapter that charges the built-in rechargeable battery pack. A wall or desk mount display bracket and mounting hardware is also included.



The VET-330WH's one-inch/25-mm high LCD display provides easy-to-read weights clearly visible up to 25 ft away. The HOLD key retains the animal's weight during movement.



Two integral wheels ease transportation of the scale to multiple weighing locations.



The scale's roomy 36 in W x 22 in D x 1.97 in H / 90 cm W x 55 cm D x 5 cm H stainless steel platform allows it to be used for a wide range of purposes.



Assembly is easy! Simply screw on the four non-skid leveling feet, install the mounting bracket to the display, and it's ready to weigh right out of the box.



The stainless steel guide handle allows you to manuever the scale wherever needed.



The scale is completely portable and can be raised vertically and rolled wherever needed via its carrying handle and wheels.



## **Versatile Digital Platform Scale**

Model	VET-330WH
Capacity	330 lb x 0.2 lb / 150 kg x 100 g
Platform Dimensions	36 in W x 22 in D x 1.97 in H / 90 cm W x 55 cm D x 5 cm H
Indicator Dimensions	10.6 in W x 3.5 in D x 1.8 in H / 270 mm W x 90 mm D x 45 mm H
Weight Display	Seven-digit with decimal, seven-segment, 1.0 inch / 25 mm LCD
Cable Length	7.0 ft / 2.1 m cable from scale base to display
Power Requirements	Built-in rechargeable battery pack 7.2V DC/1200mA (DETECTO part no. 6800- 1067) or the included 120 to 230V AC 15V DC 300 mA multi-Blade input wall plug-in UL/CUL listed AC power adapter (DETECTO part no. 6800-1044)
Construction	Stainless steel
Tare	100% of scale capacity
Keys	HOLD, UNIT, ON/ZERO/OFF, TARE
Operating Environment Temperature Range / Humidity	40 to 105 °F (+5 to +35 °C) / 25% ~95% RH
Net Weight	40 lb / 18 kg
Shipping Weight	50 lb / 23 kg
UPC Code	809161146301

### **VET-330WH Included Items:**



DETECTO Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

----



 102 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001 www.Detecto.com

## **SOLD BY:**

© Copyright 2022 Cardinal Scale Mfg.Co. • Printed in USA • CAR/00/0522/C220A