

VET330WH

DIGITAL VETERINARY SCALE

**Ideal for
Veterinarians,
Pet Stores,
Kennels or
Laboratories.**

- Integral Wheels/Carrying Handle
- Rechargeable Battery Power
- 330 Pound Capacity



Detecto's economical, portable VET330WH digital scale with remote weight display is perfect for weighing medium to large size animals.





VET330WH

DIGITAL VETERINARY SCALE

Featuring Built-in Wheels
and Carrying Handle



Detecto's VET330WH features rechargeable battery power and comes with an AC adapter. The scale's easy-to-read remote LCD weight display features 1-inch-high digits. The VET330WH features a low battery indicator and power-saving selectable timing auto shut-off. The stainless steel enclosure provides long-lasting durability and easy



wipe-down cleaning. The HOLD key on the remote display allows the animal's weight to be retained after leaving the platform. The remote display's lengthy scale cord allows it to be positioned away from the platform or nearby, wherever weight display is convenient for viewing.



MOBILITY AND VERSATILITY!

The VET330WH's rechargeable battery power combined with its light weight and integral wheels/ carrying handle allow the scale to be used in numerous locations, and the weighing can come to the patient.

VET330WH SPECIFICATIONS

CAPACITY:	330 lb x 0.2 lb / 150 kg x 0.1 kg
PLATFORM DIMENSIONS:	36" W x 22" D x 1.97" H (91 cm W x 56 cm D x 5 cm H)
DISPLAY DIMENSIONS:	10.63" W x 3.54" D x 1.77" H (27 cm W x 9 cm D x 4 cm H)
WEIGHT DISPLAY:	LCD, 7 Digit, with decimal 1.0" characters
SHIPPING WEIGHT:	50 lb / 23 kg
POWER:	Built-in rechargeable battery pack (7.2VDC/1200mA) or 120 VAC 60Hz powering a 15 VDC 300 mA wall plug-in UL/CUL listed power supply
CONSTRUCTION:	Stainless steel
TARE:	100% of full scale capacity
KEYS:	HOLD, UNIT, ON/ZERO/OFF, TARE
TEMPERATURE:	40° to 105° F (5° to 35° C)
HUMIDITY:	25% ~ 95% RH

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

 **DETECTO**[®]
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Sold By: