



EYE-LEVEL WEIGHBEAM

PHYSICIAN SCALES



**DETECTO is the name brand that means
Scales in the healthcare industry since 1900.**

- Trusted by physicians worldwide since 1900
- Mild steel construction with powder coat paint for long-lasting quality
- Available in lb, lb/kg, and kg models
- Height rod, wheels, and handpost combinations available

MEDICAL WEIGHING

PRECISION SINCE 1900



DETECTO is the name brand that means Scales in the healthcare industry since 1900. DETECTO's high-quality weighbeam medical scale is globally recognized as the leader in the mechanical physician scale market.

The innovatively-designed scale features readings from either side of the weighbeam and is imprinted on both sides, allowing the nurse and patient to both view weight simultaneously. Pounds, kilograms, or pounds/kilograms are all available on the easy-to-read beam.

The dual-reading height rod measures in inches and centimeters at a glance (24-84 inch / 60-212 cm). The slip-resistant black platform cover may be removed for easy cleaning.



When the patient steps on the scale, simply move the large poise and then the small poise until the weighbeam balances. It's that simple!

- 13 different models available
- Capacity of 440 lb/200 kg, 450 lb, or 200 kg
- Simple to operate – no power required
- Easy to read due to eye-level display
- Sturdy construction with mild steel frame and powder coat paint finish
- Height measurement range: 24-84 inch / 60-212 cm
- Height graduation: 0.125 inch / 0.1 cm



The slip-resistant steel base and anti-tip column design allows patients to be safely measured and weighed.

SPECIFY "DETECTO"

IT'S THE NAME THAT MEANS SCALES IN THE HEALTHCARE INDUSTRY.



The die-cast weighbeam is easy to read from both sides and makes weighing fast and simple.

Both the nurse and patient can view the weight at the same time since it is imprinted on both sides.

DETECTO's eye-level weighbeam scales can be used in hospitals, clinics, doctor's offices, gyms, and even home use.

ACCESSORIES

DETECTO model 448 physician scale with handpost, height rod, and wheels



EASY MOBILITY

Wheels are included on some models (see page 4) for easy mobility when changing locations.



HEIGHT ROD

The height rod is easy to extend and measures in inches and centimeters: 24-84 inch / 60-212 cm



HANDPOST

An important aid for elderly or unsteady patients. The handpost is bolted to the platform, so the patient can steady themselves on it while weighing.



WHEELS

Two rear wheels on some models allow easy transport. The light weight of the DETECTO physician scales combined with wheels makes transport simple.



PLATFORM COVER

Slip-resistant platform cover that may be removed for easy cleaning



EYE-LEVEL WEIGHBEAM PHYSICIAN SCALES



Model	437	439	438	449	448	2371	2391	2381	2491
Capacity	450 lb x 4 oz					200 kg x 100 g			
Height Rod		•	•	•	•		•	•	•
Wheels			•		•			•	
Handpost				•	•				•
Shipping Weight	35 lb / 15.8 kg	38 lb / 17.2 kg	39 lb / 17.6 kg	43 lb / 19.5 kg	43 lb / 19.5 kg	35 lb / 15.8 kg	37 lb / 16.7 kg	39 lb / 17.6 kg	43 lb / 19.5 kg
UPC Code	809161133301	809161133103	809161133509	809161133806	809161133707	809161133400	809161133202	809161133608	809161133905



Model	337	339	338	349
Capacity	440 lb x 4 oz or 200 kg x 100 g			
Height Rod		•	•	•
Wheels			•	
Handpost				•
Shipping Weight	35 lb / 15.8 kg	38 lb / 17.2 kg	39 lb / 17.6 kg	43 lb / 19.5 kg
UPC Code	809161134100	809161134001	809161134209	809161134308

Height Measurement Range: 24-84 in (60 - 212 cm)

Height Graduation: 0.125 in (0.1 cm)

Platform Size: 14.75 in x 10.75 in (37.5 x 27.3 cm)

Overall Width: 17.5 in / 44.5 cm

Overall Height: 58.25 in / 148 cm

Overall Depth (with wheels): 23.25 in / 59 cm

Construction: Mild steel with powder coat paint finish

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.



102 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

SOLD BY: