

DETE

CTO



# DETECTO'S SOLO® SCALE STANDS ALONE IN ECONOMIC CLINICAL WEIGHING

- 550 lb x 0.2 lb / 250 kg x 0.1 kg Capacity
- Weight, Height, and Body Mass Index
- 14 in W x 15 in D x 2.5 in H Platform Size
- 6 AA Batteries (included) or Optional AC Adapter
- Inline Mechanical Height Rod
- USB Port for EMR/EHR Connectivity
- Mother/Baby Mode Factors the Weight of Infants Held by an Adult



# UNMATCHED ECONOMICAL WEIGHING QUALITY



Eye-level height display for ease of use (52 in / 132 cm H)

Inline height rod measures: 30 in to 84.5 in / 76 cm to 214.5 cm

4000

1-inch/25-mm high LCD, clearly-visible display



OneWeigh units locking into Pounds or Kilograms avoids dosage errors based on incorrect units.

Rock-solid, stable column – aluminum white powder coated

The Name Brand Preferred by Clinical Professionals Around the World

Two wheels for quick and easy portable transport

Rubber mat for patient assurance and stability





## **Stylish Indicator**

The solo<sup>®</sup> indicator features 8 simple buttons that are easily identifiable in any language.



### **EMR/EHR-Ready**

The scale's USB port allows you to output measurement data to EMR/EHR software.



# **Inline Height Rod**

The mechanical height rod is situated to the side of the LCD, so it never interferes with viewing measurements.



# **AC Adapter**

The AC port is located in the base, so no pesky wires dangle from up on the indicator (AC Adapter optional: model PD-AC).



#### **Battery Power**

No power outlet available, no problem! The solo<sup>®</sup> comes with 6 AA batteries for weighing anywhere.



# **Body Mass Index**

Simply step on the scale to weigh and the scroll arrows let you rapidly enter your height to factor BMI.



# Wide, Flat Platform

The oversized 14 x 15 inch platform with rubber mat is only 2.4 inches high making it easy for patients to step on/off.



Portable, Wheels

Two integral wheels allow for ease of movement whenever you need to transport the solo<sup>®</sup>.



DIGITAL CLINICAL SCALES	
Model	solo
Capacity	550 lb x 0.2 lb / 250 kg x 0.1 kg
Platform Size	14 in W x 15 in D x 2.4 in H / 36 cm W x 38 cm D x 6 cm H
Weight Units	Pounds or Kilograms (selectable upon initialization at start-up; OneWeigh units locking to avoid dosage errors)
Power	Six (6) AA size Alkaline batteries (included) OR an optional medical device 9V AC/DC wall plug-in adapter (Detecto part number PD-AC). The AC adapter comes with model SOLO-AC.
Display	Five-digit, seven-segment, 1.0-inch (25-mm) high LCD digits
Operation Temperature	50 to 104 °F / +10 to +40 °C
	ON/OFF, ZERO, BMI HEIGHT, ENTER, LOCK/RELEASE, PRINT, UP ARROW/Mother and Baby, DOWN ARROW Interface: USB (ASCII Protocol)
Annunciators	ZERO, LOCK, STA, HOLD, LB, KG, CM, FT/IN, HEIGHT, WEIGHT, BMI, AM, PM
Scale Dimensions	14.3 in W x 19.5 in D x 52.0 in H / 36.3 cm W x 49.5 cm D x 132.1 cm H
Interface	USB (ASCII Protocol)
Height Rod	Mechanical included standard, inline
Height Rod Units	Inches or Centimeters
Height Rod Range	30 in to 84.5 in / 76 cm to 214.5 cm
Manual Input Height	2 ft 0.0 in – 7 ft 10.5 in by 0.5 in / 60.0 cm - 240.0 cm by 1cm
Minimum Weighing	11 lb / 5 kg
Wheels	(2) wheels for ease of movement
Selectable Auto Shut-Off	Up to 9 Minutes
Construction	Column and height rod – aluminum alloy white powder coated; indicator enclosure, height rod arm, platform and base – white ABS; mat – black PVC
Product Weight	17.7 lb / 7.7 kg
Shipping Weight	25 lb / 11.3 kg
Package Dimen- sions	
UPC Code	809161197709

DETECTO's SOLO-RI is also available with an eye-level column as model SOLO

1000

DHTHCTO

Model SOLO

0.0

DET

E C T

0

Model SOLO-RI



102 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153 www.Detecto.com **SOLD BY:**