

# solo®

DIGITAL CLINICAL SCALES

## DIGITAL PHYSICIAN SCALE WITH HEIGHT ROD




- 550 lb x 0.2 lb / 250 kg x 0.1 kg Capacity
- Weight, Height, and Body Mass Index
- 14 in W x 15 in D x 2.5 in H Platform Size
- 6 AA Batteries (included) or Optional AC Adapter
- Inline Mechanical Height Rod

**THE MOST ECONOMICAL,  
HIGH-QUALITY DIGITAL  
PHYSICIAN SCALE  
ON THE MARKET!**

*Detecto's solo® scale stands  
alone in economic clinical  
weighing*

The oversized 14 x 15  
inch platform with rubber  
mat is only 2.4 inches  
high making it easy for  
patients to step on/off.

*No power outlet available, no  
problem! The solo® comes with 6 AA  
batteries for weighing anywhere.*



The solo® indicator  
features 8 simple buttons  
that are easily identifiable  
in any language.



*The AC port is  
located in the base,  
so no pesky wires  
dangle from up on  
the indicator (AC Adapter  
optional: model PD-AC)*



# DETECTO®

[www.Detecto.com](http://www.Detecto.com)